

PLANTAR FASCIITIS

AVERAGE DAYS WEAR



APPLICATION OVERVIEW

The plantar fascia, or arch tendon, is a thick band of connective tissue that runs from the heel to the front of the foot. Arch pain may be triggered by overtraining, poor-fitting or worn shoes, weight gain, over pronation, or abnormal foot structure. KT Tape helps treat this condition by offering support and resting the tissues.

WHAT YOU NEED

3 strips of KT TAPE

BEFORE YOU START

APPLY BEFORE ACTIVITY

Apply one hour before beginning activity

CLEAN SKIN

Clean dirt, oils and lotions from area

ACTIVATE ADHESIVE


After application rub tape vigorously to activate adhesive



BODY POSITION

Point toes up toward shin.

STRIP ONE



0% STRETCH

ANCHOR: Anchor full strip on ball of foot as shown without stretch.



50% STRETCH

APPLY: Apply strip along bottom of foot and up the back of heel as shown with 50% stretch on tape.



0% STRETCH


FINISH: Lay the end down without stretch. Smooth tape against the arch.

STRIP TWO




0% STRETCH

ANCHOR: Anchor a second full strip four inches above the inside of the ankle without stretch.



80% STRETCH

APPLY: Apply tape around bottom of heel and up the outside of ankle with 80% stretch on tape.



0% STRETCH

FINISH: Apply last two inches of tape without stretch.

STRIP THREE




0% STRETCH

ANCHOR: Anchor a full strip slightly behind the first strip angled toward the arch without stretch.



80% STRETCH

APPLY: Apply tape across arch with 80% stretch.



0% STRETCH

FINISH: Apply last two inches of tape on outside of foot without stretch.



WATCH THE VIDEO

kttape.com/instructions/plantar-fasciitis